



Technical Bulletin

Ref: **TB- 08004**

Issue date: **March 2010**

Product: INOmax[®] DS / INOblender[®]

Priority: Low

Affected parts: INOmax[®] DS / INOblender[®]

Classification: Information

Subject: **Use of Cleaning Agents**

Please read the INOmax DS and INOblender manual for a complete list of cleaning procedures, cautions and warnings.

Caution: Apply cleaning agent to a cloth before application; do not spray directly on the device. It is important to prevent pooling and direct contact with electrical connections which can cause damage over time.

Clean the outer surface of the INOmax DS or INOblender with a soft cloth dampened in a mild soap and water solution, isopropyl alcohol (70%) or with one of the following cleaning agents while following the manufacturer's recommendations.

Cleaning Agent	Active Ingredients
Precise Hospital Foam Cleaner Disinfectant by Caltech Industries	o-Phenylphenol < 0.37% Other ingredients 99.63%
Pure Green 24 by Pure Green, LLC	SDC – silver ions 0.003% Citric acid 4.84% Other ingredients 95.157%
PDI Super Sani Cloth by PDI	n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.25% n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.25% Isopropyl alcohol 55% Inert ingredients 44.50%
Sani Cloth HB by PDI	n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.07% n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.07% Inert ingredients 99.86%
Asepti-HB by Ecolab Inc.	n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.07% n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.07% Inert ingredients 99.86%
Cavicide and CaviWipes by Metrex	Dilsobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride 0.28% Isopropyl alcohol 17.2% Inert ingredients 82.52%
EMPOWER- dual enzymatic detergent by Metrex	Proteinase subtilisin < 2% Propylene glycol water

If you have further questions, please contact technical support at 1-877-566-9466.