



Technical Bulletin

Ref: **TB- 08003**

Issue date: **April 2008**

Product: INOmax™ DS

Priority:

Affected parts: Non-Ventilator Applications

Classification: Information

Subject: **INOmax DS testing standards**

The INOmax DS has been tested to the following standards:

| Test | Standards |
|------------------------|---|
| Electrical safety | UL 60601-1, EN 60601-1, IEC 60601-1, and CSA (Canadian) STD C22.2 |
| Emissions and Immunity | IEC 60601-1-2:2001 IEC 61000-3-2:2000 IEC 61000-3-3:1995 EN 60601-1-2:2001 EN 61000-3-2:2000 EN 61000-3-3:1995 |
| Shock | FDA 7.1.4 |
| Sinusoidal Vibration | FDA 7.1.4 |
| Random Vibration | FDA 7.1.4 |
| Shock | CEN/TS 14507-1:2003 (E) |
| Sinusoidal Vibration | CEN/TS 14507-1:2003 (E) |
| Random Vibration | CEN/TS 14507-1:2003 (E) |
| Bump Test | CEN/TS 14507-1:2003 (E) |
| Water Egress | IPX1 |

If you have further questions, please contact technical support at 1-877-566-9466.